

SPC Hip Stem

CATALOGUE



CONTENTS

IMPLANTS	3
SPC Hip Stem CCD 125° (ISO 5832-11 Ti6Al7Nb)	3
SPC Hip Stem CCD 135° (ISO 5832-11 Ti6Al7Nb)	4
SPC Hip Stem CCD 145° (ISO 5832-11 Ti6Al7Nb)	4
Femoral Heads	5
INSTRUMENTS	6
Basic Instrument Set, modular (Art. no. 367-116)	6
SPC Hip Stem Rasp Instrument Set (Art. no. 271-324)	7
X-ray templates	8
Optional accessories	8

IMPLANTS




Figure 1: SPC Hip Stem

SPC Hip Stem CCD 125° (ISO 5832-11 Ti6Al7Nb)

Implant		Length [mm]	CCD angle	Offset [mm]	Art. no.
SPC Hip Stem Size 5	ISO 5832-11 Ti6Al7Nb, cementless	135.6	125°	37.5	271-205
SPC Hip Stem Size 6	ISO 5832-11 Ti6Al7Nb, cementless	139.2	125°	38.8	271-206
SPC Hip Stem Size 7	ISO 5832-11 Ti6Al7Nb, cementless	142.8	125°	40.1	271-207
SPC Hip Stem Size 8	ISO 5832-11 Ti6Al7Nb, cementless	146.4	125°	41.5	271-208
SPC Hip Stem Size 9	ISO 5832-11 Ti6Al7Nb, cementless	150	125°	42.8	271-209
SPC Hip Stem Size 10	ISO 5832-11 Ti6Al7Nb, cementless	153.6	125°	44.1	271-210
SPC Hip Stem Size 11.25	ISO 5832-11 Ti6Al7Nb, cementless	158.1	125°	45.5	271-211
SPC Hip Stem Size 12.5	ISO 5832-11 Ti6Al7Nb, cementless	162.6	125°	47	271-212
SPC Hip Stem Size 13.75	ISO 5832-11 Ti6Al7Nb, cementless	167.1	125°	48.4	271-213
SPC Hip Stem Size 15	ISO 5832-11 Ti6Al7Nb, cementless	171.6	125°	49.9	271-215
SPC Hip Stem Size 16.25	ISO 5832-11 Ti6Al7Nb, cementless	176.1	125°	51.3	271-216
SPC Hip Stem Size 17.5	ISO 5832-11 Ti6Al7Nb, cementless	180.6	125°	52.8	271-217
SPC Hip Stem Size 20	ISO 5832-11 Ti6Al7Nb, cementless	189.6	125°	55.7	271-220

SPC Hip Stem CCD 135° (ISO 5832-11 Ti6Al7Nb)

Implant		Length [mm]	CCD angle	Offset [mm]	Art. no.
SPC Hip Stem Size 5	ISO 5832-11 Ti6Al7Nb, cementless	135.6	135°	35.1	271-105
SPC Hip Stem Size 6	ISO 5832-11 Ti6Al7Nb, cementless	139.2	135°	36.3	271-106
SPC Hip Stem Size 7	ISO 5832-11 Ti6Al7Nb, cementless	142.8	135°	37.6	271-107
SPC Hip Stem Size 8	ISO 5832-11 Ti6Al7Nb, cementless	146.4	135°	38.8	271-108
SPC Hip Stem Size 9	ISO 5832-11 Ti6Al7Nb, cementless	150	135°	40.1	271-109
SPC Hip Stem Size 10	ISO 5832-11 Ti6Al7Nb, cementless	153.6	135°	41.2	271-110
SPC Hip Stem Size 11.25	ISO 5832-11 Ti6Al7Nb, cementless	158.1	135°	42.6	271-111
SPC Hip Stem Size 12.5	ISO 5832-11 Ti6Al7Nb, cementless	162.6	135°	43.9	271-112
SPC Hip Stem Size 13.75	ISO 5832-11 Ti6Al7Nb, cementless	167.1	135°	45.3	271-113
SPC Hip Stem Size 15	ISO 5832-11 Ti6Al7Nb, cementless	171.6	135°	46.6	271-115
SPC Hip Stem Size 16.25	ISO 5832-11 Ti6Al7Nb, cementless	176.1	135°	47.9	271-116
SPC Hip Stem Size 17.5	ISO 5832-11 Ti6Al7Nb, cementless	180.6	135°	49.2	271-117
SPC Hip Stem Size 20	ISO 5832-11 Ti6Al7Nb, cementless	189.6	135°	51.9	271-120

SPC Hip Stem CCD 145° (ISO 5832-11 Ti6Al7Nb)

Implant		Length [mm]	CCD angle	Offset [mm]	Art. no.
SPC Hip Stem Size 5	ISO 5832-11 Ti6Al7Nb, cementless	135.6	145°	32.8	271-005
SPC Hip Stem Size 6	ISO 5832-11 Ti6Al7Nb, cementless	139.2	145°	33.9	271-006
SPC Hip Stem Size 7	ISO 5832-11 Ti6Al7Nb, cementless	142.8	145°	35	271-007
SPC Hip Stem Size 8	ISO 5832-11 Ti6Al7Nb, cementless	146.4	145°	36.1	271-008
SPC Hip Stem Size 9	ISO 5832-11 Ti6Al7Nb, cementless	150	145°	37.2	271-009
SPC Hip Stem Size 10	ISO 5832-11 Ti6Al7Nb, cementless	153.6	145°	38.2	271-010
SPC Hip Stem Size 11.25	ISO 5832-11 Ti6Al7Nb, cementless	158.1	145°	39.4	271-011
SPC Hip Stem Size 12.5	ISO 5832-11 Ti6Al7Nb, cementless	162.6	145°	40.6	271-012
SPC Hip Stem Size 13.75	ISO 5832-11 Ti6Al7Nb, cementless	167.1	145°	41.8	271-013
SPC Hip Stem Size 15	ISO 5832-11 Ti6Al7Nb, cementless	171.6	145°	43	271-015
SPC Hip Stem Size 16.25	ISO 5832-11 Ti6Al7Nb, cementless	176.1	145°	44.2	271-016
SPC Hip Stem Size 17.5	ISO 5832-11 Ti6Al7Nb, cementless	180.6	145°	45.4	271-017
SPC Hip Stem Size 20	ISO 5832-11 Ti6Al7Nb, cementless	189.6	145°	47.8	271-020

Femoral Heads




Figure 2: Femoral Heads

Implant	Nominal diameter	Neck length	Art. no.
Femoral Head 12/14 ISO 5832-12 CoCrMo	Ø 22 mm	S (-4 mm) M (0 mm) L (+4 mm)	030-2200 to 030-2202
Femoral Head 12/14 ISO 5832-9 implant steel	Ø 28 mm	S (-3.5 mm) M (0 mm) L (+3.5 mm)	020-2800 to 020-2804
Femoral Head 12/14 ISO 5832-12 CoCrMo		XL (+7 mm) XXL (+10.5 mm)	030-2800 to 030-2804
BioloX® delta Femoral Head 12/14 ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - Ceramic	Ø 28 mm	S (-3.5 mm) M (0 mm) L (+3.5 mm)	367-1140 to 367-1142
BioloX® forte Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic			367-907 to 367- 909
ELEC® Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic			384-001 to 384- 003
Femoral Head 12/14 ISO 5832-9 implant steel	Ø 32 mm	S (-4 mm) M (0 mm) L (+4 mm)	020-3200 to 020-3204
Femoral Head 12/14 ISO 5832-12 CoCrMo		XL (+8 mm) XXL (+12 mm)	030-3200 to 030-3204
BioloX® delta Femoral Head 12/14 ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - Ceramic	Ø 32 mm	S (-4 mm) M (0 mm) L (+4 mm) XL (+8 mm)	367-1143 to 367-1145 and 367-1149
BioloX® forte Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic	Ø 32 mm	S (-4 mm) M (0 mm) L (+4 mm)	367-910 to 367- 912
ELEC® Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic			384-004 to 384- 006
BioloX® delta Femoral Head 12/14 ISO 6474-2 Al ₂ O ₃ / ZrO ₂ - Ceramic	Ø 36 mm	S (-4 mm) M (0 mm) L (+4 mm) XL (+8 mm)	367-1146 to 367-1148 and 367-1150
BioloX® forte Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic	Ø 36 mm	S (-4 mm) M (0 mm) L (+4 mm)	367-930 to 367- 932
ELEC® Femoral Head 12/14 ISO 6474-1 Al ₂ O ₃ - Ceramic			384-007 to 384- 009

INSTRUMENTS

Basic Instrument Set, modular (Art. no. 367-116)




Figure 3: Basic Instrument Set, modular

No.	Instrument	Art. no.
1	SPC-Stem Impactor with Silicone Handle grey, L=300mm	271-301
2	Rasp Handle Straight L=253mm	271-321
3	Stem Extractor	271-322
4	Guide Rod Ø8mm L=120mm	506-015
5	Femoral Rasp with Silicone Handle grey, L = 224mm	506-071
6a	Hammer 1100g with Silicone Handle grey, L=250mm	506-072
6b	Impact Head for Hammer with Silicone Handle	506-072 II/II
7	Box Chisel Small 6mm x 20mm	506-091
8	Modular Weight	506-2107
9	Tray Basic Instrument Set modular	367-186

SPC Hip Stem Rasp Instrument Set (Art. no. 271-324)




Figure 4: SPC Hip Stem Rasp Instrument Set

No.	Instrument	Art. no.
1	Femoral Head Impactor Handle with Silicone Handle grey, L=192mm	506-060
2	Head End Ø28/32mm for Femoral Head Impactor Handle	506-062
3	Trial Head Ø28 12/14 S POM grey	512-280
4	Trial Head Ø28 12/14 M POM green	512-281
5	Trial Head Ø28 12/14 L POM blue	512-282
6	Trial Head Ø28 12/14 XL POM red	512-283
7	Trial Head Ø28 12/14 XXL POM brown	512-284
8	Trial Head Ø32 12/14 S POM grey	512-320
9	Trial Head Ø32 12/14 M POM green	512-321
10	Trial Head Ø32 12/14 L POM blue	512-322
11	Trial Head Ø32 12/14 XL POM red	512-323
12	Trial Head Ø32 12/14 XXL POM brown	512-324
13	Trial Cone SPC 125°	271-302
14	Trial Cone SPC 135°	271-303
15	Trial Cone SPC 145°	271-304
16	Rasp Size 5 SPC	271-305
17	Rasp Size 6 SPC	271-306
18	Rasp Size 7 SPC	271-307

No.	Instrument	Art. no.
19	Rasp Size 8 SPC	271-308
20	Rasp Size 9 SPC	271-309
21	Rasp Size 10 SPC	271-310
22	Rasp Size 11,25 SPC	271-311
23	Rasp Size 12,5 SPC	271-312
24	Rasp Size 13,75 SPC	271-313
25	Rasp Size 15 SPC	271-315
26	Rasp Size 16,25 SPC	271-316
27	Rasp Size 17,5 SPC	271-317
28	Rasp Size 20 SPC	271-320
29	Tray SPC Rasp Instruments	271-323

X-ray templates

Designation	Art. no.
X-ray templates SPC CCD 125°	271-125
X-ray templates SPC CCD 135°	271-135
X-ray templates SPC CCD 145°	271-145

Optional accessories

Designation	Art. no.
Femoral Neck Resection Guide CCD 135°	367-006
Head End Ø36mm for Femoral Head Impactor Handle	506-1212
Trial Head Ø36 12/14 S POM grey	512-361
Trial Head Ø36 12/14 M POM green	512-362
Trial Head Ø36 12/14 L POM blue	512-363
Trial Head Ø36 12/14 XL POM red	512-364



OHST Medizintechnik AG

Grünauer Fenn 3

14712 Rathenow

Germany

Fon +49 (0) 3385 5420 0

Fax +49 (0) 3385 5420 99

E-Mail info@ohst.de

URL www.ohst.de